



ENVIRONMENTAL HEALTH & SAFETY

SAFETY GUIDELINE Autoclaves

This document is intended to help you prevent the occurrence of injury or other incident during autoclave use. **Autoclaves** use high heat and pressurized steam to destroy microorganisms.

Burn Hazards!

High Temperature:

- Steam
- Surfaces
- Items
- Liquids



Explosion Hazards!

Pressure and Heat plus:

- Air Tight Containers
- Improper Loading
- Door Malfunctions
- Oxidizers



Other Hazards!

Physical Hazards:

- Heavy Lifting
 - Pinch, Impact, Cut
- Inhalation or Contact:
- Chemicals
 - Biologicals

Crucial Considerations

- *Training should be provided* before any new user is allowed to operate the autoclave on their own. Always document training events.
- *Strictly follow* the manufacturer's *manual/instructions* along with departmental, program, or laboratory *specific procedures* associated with the autoclave itself and the materials you wish to autoclave.



- *Always wear personal protective equipment (PPE)* appropriate for the hazards. Examples include: goggles or face shields, nitrile gloves, lab coats, *thermal mitts* (pictured here), and closed-toe shoes.
- *Report all error messages or malfunctions* to your supervisor and whomever oversees the autoclave.

Post an "Out-of-Service" sign to prevent use until repairs are made.

- *Ensure materials to be autoclaved are compatible* with peak autoclave temperatures and pressures. Polyethylene is not compatible with autoclaving, as seen in this image of a melted five gallon carboy. The table below lists examples of non-compatible materials. Unsure if your material is compatible? Simply contact UNMC EHS at 402-559-6356 or unmcchs@unmc.edu for additional guidance.



Examples of Materials that are Non-Compatible with Autoclaving			
Organic Solvents	Flammables	Corrosives	Sulphates
Toxic, Radioactive or Volatile Materials		Non-Stainless Steel	Nylon
Polyurethane	Polyethylene <i>(pictured above)</i>	Hypochlorite (Bleach), Chlorides, Chlorine	
Polystyrene	Polyvinyl Chloride	Acrylic	Dry Combustibles
Carcinogens or Mutagens		Strong Acids or Bases	
<i>Note: Some biohazard bags, laboratory bins, and other lab plastics are not autoclavable. Dry combustible goods (gloves, paper, etc.) may burn on dry autoclave cycles in the absence of liquids.</i>			

Loading and Unloading Best Practices

Before Loading or Unloading

- Check status indicators and messages.
- Confirm chamber pressure is zero.

While Loading or Unloading

- Utilize PPE. *Stand away from the door and open it slowly.* Steam, heat, and potentially hazardous fumes may rush out.

Content of Autoclave Loads

- Use *secondary containment* to avoid compromising the autoclave interior due to leaks, spills, or melting.
- Liquid containers should not be filled beyond $\approx 50\%$ of their total volume, otherwise they are at risk for boiling over (pictured on right).
- Use lead free autoclave indicator tape [Hazardous Material Fact Sheet, Autoclave Tape PDF](#)



- *Pressure build up within tightly sealed items can cause explosion or implosion* (pictured on left). Caps/lids or closures should be secure but also loose enough to allow air flow.
- *Avoid stacking and overloading/overpacking.* Always leave space between items and autoclave walls. Steam must be able to penetrate all areas of the autoclave load to properly sterilize it entirely.
- *Promptly unload your materials* at the end of the autoclave process to free up the equipment for others.
- *Avoid unloading autoclaves if the contents within are of unknown identity to you.* Unknowns pose a risk for potential injury and some materials may require special handling or additional processing immediately post autoclaving.

Questions? More Information? Contact: UNMC EHS at 402-559-6356 or unmcehs@unmc.edu