


## Instructions to use Vector NTI with new dynamic license.

1. If you don't have vector NTI installed, download it from <http://www.lifetechnologies.com/us/en/home/life-science/cloning/vector-nti-software/vector-nti-express-software/download-vector-nti-express.html> and follow the installation instruction.
2. If you already have Vector NTI installed, follow the following steps to activate the license.
  - a. Open License manager. Open VNTI -> Help -> License Manager.



- b. Select "User Information" tab and enter your information.


The screenshot shows the License Manager dialog box with the "User Information" tab selected. The dialog box contains the following fields:

- Name: Sanjit Pandey
- Organization: University of Nebraska Medical Center
- Phone number: 4025595925
- E-mail address: sanjit.pandey@unmc.edu

Please note:  
All fields are required. The name provided will be used as your Author name for records you create in the Local Database. An email address may be required to utilize some online bioinformatics resources.

Close

- c. Select “Applications” tab in the license manager window. Make sure Vector NTI Express is set to “Dynamic License”.



- d. Click on “Dynamic” button to get to the Dynamic License window. Fill the form with your information.

Dynamic License

Name: Sanjit Pandey

Organizations: University of Nebraska Medical Center

Phone number: 4025595925

E-mail address: sanjit.pandey@unmc.edu

URL of DLS: Send an email to request this

DLS server requires authentication

Login:

Password:

Internet Connection Settings... Test Connection...

Apply Cancel

- e. For URL of DLS send an email to sanjit.pandey@unmc.edu.
- f. Click "Test Connection" button and click on the "Connect" button in the following window. You should see the "Connect successful" message on the right.

Dynamic Licenses server connection tester

License Server URL (Ask your administrator for details)

Send an email to request the url

Use DLS authorization (if applicable)

User name:

Password:


Internet Connection Settings...

Connection diagnostics

Connect successful.  
VectorNTI Dynamic License Server version 2.0 (Aug 27)  
Available :  
Vector NTI Express

Connect Clear Close

- g. Close the window and click “Apply” in the Dynamic license window. When successfully registered, you should get the following message.



- h. You should be good to go once the application restarts.